

3 ways to transfer from Sonocent Recorder to Audio Notetaker



Option 1: WiFi



Press **'Transfer'** on the Recorder app.

In Audio Notetaker go to **Menu -> manage -> portable devices**. Click your device and enter 4-digit transfer code displayed in app.

Option 2: USB



Connect your mobile device to your computer via **USB**.

- **Android:** Open the device in **File Explorer** and copy your files.
- **iOS:** Select your device in iTunes and use **Apps, File Sharing** to transfer your files

Option 3: Online



Select a file in Recorder and press **'share'**. Choose sharing method (i.e. Dropbox, Google Drive, OneDrive, iCloud, AirDrop).

Retrieve file from online hosting.

